



### SIZE OPTIONS

The N.hance® Access sling is available in a range of standard sizes from Extra Small to Extra Extra Large.



### EASY ATTACHMENT

The N.hance® Access sling has a choice of multiple coloured loops for different positions and easy identification.



### COMFORT & SUPPORT

Generous padded underarm rolls and a non slip belt make the N.hance® Access sling very comfortable and supportive.

## N.HANCE® ACCESS SLING – MAKING TOILETING SIMPLE

### A comfortable yet supportive sling for toileting and clothing removal

With generously padded underarm rolls and padded legs, the N.hance® Access sling makes removing lower body clothing for toileting whilst being hoisted straightforward and safe.

The large anti-slip panel in the back of the sling helps users maintain an upright posture in the sling whilst the clip and buckle belt maintains position and provides security for the user.

The N.hance® Access sling is available in a range of sizes from Small to Extra Extra Large. Users must have good upper body strength to use the N.hance® Access sling, and be able to sit unaided.

**COMFORTABLE**

**SUPPORTIVE**

**VERSATILE**




# N.HANCE® ACCESS TOILETING SLING



**WARRANTY**  
1 year



**IN STOCK**  
Yes



**SWL**  
200kg



Access Sling

## PRODUCT INFO


Code	Description
NH2012S	N.hance® Access Toileting Sling - Small
NH2012M	N.hance® Access Toileting Sling - Medium
NH2012L	N.hance® Access Toileting Sling - Large
NH2012XL	N.hance® Access Toileting Sling - Extra Large
NH2012XXL	N.hance® Access Toileting Sling - Extra Extra Large

Highlighted rows are special order items

**Disclaimer:** Prior to use, a documented risk assessment and care plan must be completed to ensure the safety of the user and carers. If you are in any doubt about the suitability of the sling, contact a qualified healthcare professional for an independent assessment. N.hance® slings are for use with 'loop type' spreader bars and should not be used on 'wishbone type' spreader bars. If you are in any doubt as to the compatibility of an N.hance® sling with the hoist that you are using, contact a qualified healthcare professional.



Classic Low Sling



Classic Sling  
with Head Support



Deluxe Spacer Sling  
with Head Support



Comfort Spacer Sling  
without Head Support



Access Sling

## N.HANCE® SLING COMPARISON CHART

Condition	Classic Low Mesh	Classic Low Poly	Classic Mesh	Classic Poly	Classic Spacer	Deluxe Spacer	Comfort Spacer	Access
User has limited upper body strength	✓	✓	✓	✓	✓	✓	✓	
User has limited head control			✓	✓	✓	✓	✓	
User requires toileting	✓	✓	✓	✓	✓	✓		✓
User requires clothing to be removed in sling								✓
User requires sling to be left underneath in chair					✓	✓	✓	
Sling used for bathing	✓		✓		✓*	✓*	✓*	
User goes into extension						✓	✓	
User has limited upper body strength						✓	✓	
User has above knee amputation							✓	

\*Not recommended for prolonged periods